

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, EDR, TRF, Fast, SS, FR, ARC</i>
Plane	<i>3-PLANE</i>	<b>SCAN RANGE</b>	
Series Description	<i>3-Plane Localizer</i>	FOV	<i>48.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>8.0</i>
TE	<i>80.0</i>	Slice Spacing	<i>5.0</i>
Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>	
TR	<i>525.0</i>	Freq	<i>384</i>
Receiver Bandwidth	<i>83.33</i>	Phase	<i>160</i>
<b>IMAGE ENHANCE</b>		Freq DIR	<i>Unswap</i>
Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>
<b>GATING/TRIGGER</b>		Phase FOV	<i>1.00</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
<b>FMRI</b>		Phase Correction	<i>No</i>
PSD Trigger	<i>Internal</i>	<b>USER CVS</b>	
View Order	<i>Bottom/Up</i>	User CV1	<i>1.00</i>
# of Repetitions REST	<i>0</i>	<b>MULTI-PHASE</b>	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
<b>SAT</b>		Mask Phase	<i>0</i>
Tag Type	<i>None</i>	Mask Pause	<i>0</i>
<b>TRICKS</b>		<b>DIFFUSION</b>	
Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>
Auto Subtract	<i>0</i>	<b>CONTRAST</b>	
Auto SCIC	<i>Off</i>	Contrast Yes/No	<i>No</i>

3-Plane Localizer

3-Plane Localizer

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

Cor SSFSE BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Cor SSFSE BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, SS, Asset</i>	
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Cor SSFSE BH</i>	FOV	<i>40.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>7.0</i>	
	TE	<i>90.0</i>	Slice Spacing	<i>1.0</i>	
	Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>		
	TR	<i>2095.6</i>	Freq	<i>360</i>	
	Receiver Bandwidth	<i>62.50</i>	Phase	<i>260</i>	
	<b>IMAGE ENHANCE</b>		Freq DIR	<i>S/I</i>	
	Filter Choice	<i>B</i>	NEX	<i>1.00</i>	
	<b>GATING/TRIGGER</b>		# of Acq. Before Pause	<i>12</i>	
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>	
	<b>FMRI</b>		Auto Shim	<i>Auto</i>	
	PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>	
	View Order	<i>Bottom/Up</i>	<b>USER CVS</b>		
	# of Repetitions REST	<i>0</i>	User CV1	<i>1.00</i>	
	# of Repetitions ACTIVE	<i>0</i>	User CV2	<i>184.00</i>	
	<b>SAT</b>		User CV13	<i>1.00</i>	
	Tag Type	<i>None</i>	User CV15	<i>1.00</i>	
	<b>TRICKS</b>		<b>MULTI-PHASE</b>		
	Pause On/Off	<i>On</i>	Seperate Series	<i>0</i>	
	Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>	
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>		
		<b>DIFFUSION</b>			
		Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP512, SS, Assét</i>
Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
Series Description	<i>Ax SSFSE BH</i>	FOV	<i>38.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>7.0</i>
TE	<i>90.0</i>	Slice Spacing	<i>1.0</i>
Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>	
TR	<i>1176.1</i>	Freq	<i>360</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>260</i>
<b>IMAGE ENHANCE</b>		Freq DIR	<i>R/L</i>
Filter Choice	<i>B</i>	NEX	<i>1.00</i>
<b>GATING/TRIGGER</b>		# of Acq. Before Pause	<i>15</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>
<b>FMRI</b>		Auto Shim	<i>On</i>
PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>
View Order	<i>Bottom/Up</i>	<b>USER CVS</b>	
# of Repetitions REST	<i>0</i>	User CV1	<i>1.00</i>
# of Repetitions ACTIVE	<i>0</i>	User CV2	<i>184.00</i>
<b>SAT</b>		User CV13	<i>1.00</i>
Tag Type	<i>None</i>	User CV15	<i>1.00</i>
<b>TRICKS</b>		<b>MULTI-PHASE</b>	
Pause On/Off	<i>On</i>	Seperate Series	<i>0</i>
Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>
		<b>DIFFUSION</b>	
		Recon All Images	<i>On</i>
		<b>CONTRAST</b>	
		Contrast Yes/No	<i>No</i>

Ax SSFSE BH

Ax SSFSE BH

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

Ax SSFSE F/S BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Ax SSFSE F/S BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, SS, Asset, FR</i>	
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Ax SSFSE F/S BH</i>	FOV	<i>38.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>7.0</i>	
	TE	<i>90.0</i>	Slice Spacing	<i>1.0</i>	
	Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>		
	TR	<i>1200.9</i>	Freq	<i>360</i>	
	Receiver Bandwidth	<i>62.50</i>	Phase	<i>260</i>	
	<b>IMAGE ENHANCE</b>		Freq DIR	<i>R/L</i>	
	Filter Choice	<i>B</i>	NEX	<i>1.00</i>	
	<b>GATING/TRIGGER</b>		# of Acq. Before Pause	<i>15</i>	
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>	
	<b>FMRI</b>		Auto Shim	<i>On</i>	
	PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>	
	View Order	<i>Bottom/Up</i>	<b>USER CVS</b>		
	# of Repetitions REST	<i>0</i>	User CV1	<i>1.00</i>	
	# of Repetitions ACTIVE	<i>0</i>	User CV2	<i>184.00</i>	
	<b>SAT</b>		User CV13	<i>1.00</i>	
	Tag Type	<i>None</i>	User CV15	<i>1.00</i>	
	Fat/Water Saturation	<i>Fat Classic</i>	<b>MULTI-PHASE</b>		
	<b>TRICKS</b>		Seperate Series	<i>0</i>	
	Pause On/Off	<i>On</i>	Mask Phase	<i>0</i>	
Auto Subtract	<i>0</i>	Mask Pause	<i>0</i>		
Auto SCIC	<i>2</i>	<b>DIFFUSION</b>			
		Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

Ax FIESTA ASPIR BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Ax FIESTA ASPIR BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, EDR, Fast, ARC</i>	
	Plane	<i>AXIAL</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Ax FIESTA ASPIR BH</i>	FOV	<i>40.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>7.0</i>	
	Flip Angle	<i>50</i>	Slice Spacing	<i>1.0</i>	
	TE	<i>Minimum</i>	<b>ACQ TIMING</b>		
	Number of Echoes	<i>1</i>	Freq	<i>192</i>	
	TI	<i>210</i>	Phase	<i>300</i>	
	Receiver Bandwidth	<i>125.00</i>	Freq DIR	<i>R/L</i>	
	<b>IMAGE ENHANCE</b>		NEX	<i>1.00</i>	
	Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>	
	<b>GATING/TRIGGER</b>		Phase FOV	<i>0.80</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>	
	<b>MULTI-PHASE</b>		Phase Correction	<i>No</i>	
	Seperate Series	<i>0</i>	<b>FMRI</b>		
	Mask Phase	<i>0</i>	PSD Trigger	<i>Internal</i>	
	Mask Pause	<i>0</i>	View Order	<i>Bottom/Up</i>	
<b>DIFFUSION</b>		# of Repetitions REST	<i>0</i>		
Recon All Images	<i>On</i>	# of Repetitions ACTIVE	<i>0</i>		
<b>CONTRAST</b>		<b>SAT</b>			
Contrast Yes/No	<i>No</i>	Tag Type	<i>None</i>		
		Fat/Water Saturation	<i>Fat Special</i>		
		<b>TRICKS</b>			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		
		<b>OTHERS</b>			
		Protocol Notes	<i>Tip: - Place a shim volume over the liver to help reduce artifacts. - TI time for ASPIR (ChemSAT) can be modified by changing the Prep Time.</i>		

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

PATIENT POSITION	
Patient Entry	Feet First
Patient Position	Supine
Coil Configuration	Body
Plane	AXIAL
Series Description	Ax DWI 3:1 B-800-500-50

SCAN TIMING	
TE	Minimum
Number of Echoes	1
Number of Shots	1

IMAGE ENHANCE	
Filter Choice	None

GATING/TRIGGER	
Auto Trigger Type	Off
Resp. Trigger Window	30

FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT	
Tag Type	None
Fat/Water Saturation	Fat

TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	2

CONTRAST	
Contrast Yes/No	No

IMAGING PARAMETERS	
Imaging Mode	2D
Pulse Sequence	Spin Echo
Imaging Options	EDR, EPI, DIFF, Asset, Nav

SCAN RANGE	
FOV	40.0
Slice Thickness	7.0
Slice Spacing	1.0

ACQ TIMING	
Freq	80
Phase	128
Freq DIR	R/L
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes

USER CVS	
User CV0	1.00
User CV5	1.00
User CV17	1.00
User CV18	1.00

MULTI-PHASE	
Seperate Series	0
Mask Phase	0
Mask Pause	0

DIFFUSION	
Optimized TE	Yes
Diffusion Directions	64
Number of Diffusion Directions	1
Number of T2 Images	0
Dual Spin Echo	Off
Recon All Images	On

TRACKER	
Tracker Length	120.0
Tracker Thickness	20.0

Ax DWI 3:1 B-800-500-50


Ax DWI 3:1 B-800-500-50

---


	<p>OTHERS</p>
--	---------------

**Protocol Notes**

*Auto Tracker placement is turned on  
Tracker will be placed automatically  
Manual Placement  
Instructions:  
Navigator Tracker  
Placement Considerations  
\* Place navigator tracker on a breath hold localizer for the optimum performance\*  
- Center of tracker should be 1.5 cm below diaphragm if localizer was acquired with expiration breath hold  
- Center of tracker should be 1.5 cm above diaphragm if localizer was acquired with inspiration breath hold  
- Center of tracker should be placed on the diaphragm if localizer was acquired free breathing  
\*Please review operator manual for additional information  
Expiration-Breath Hold Tracker Placement*







*Inspiration-Breath Hold Tracker Placement*




*Free Breathing Tracker Placement*







## Easy Tro



-> No need to check anything since good waveform reflects respiratory status

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

3D Ax LAVA Flex Pre/AV	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, MPhVar, ARC, Flex</i>
	Plane	<i>AXIAL</i>	IDEAL	<i>15</i>
	Series Description	<i>3D Ax LAVA Flex Pre/AV</i>	<b>SCAN RANGE</b>	
	<b>SCAN TIMING</b>		FOV	<i>42.0</i>
	Flip Angle	<i>10</i>	Slice Thickness	<i>4.4</i>
	TE	<i>Minimum</i>	Location per Slab	<i>56</i>
	Number of Echoes	<i>2</i>	Overlap Locations	<i>0</i>
	Receiver Bandwidth	<i>142.86</i>	<b>ACQ TIMING</b>	
	<b>IMAGE ENHANCE</b>		Freq	<i>300</i>
	Filter Choice	<i>A</i>	Phase	<i>200</i>
	<b>GATING/TRIGGER</b>		Freq DIR	<i>R/L</i>
	Auto Trigger Type	<i>Off</i>	NEX	<i>1.00</i>
	<b>FMRI</b>		Phase FOV	<i>1.10</i>
	PSD Trigger	<i>Internal</i>	Auto Shim	<i>Auto</i>
	View Order	<i>Bottom/Up</i>	Phase Correction	<i>No</i>
	# of Repetitions REST	<i>0</i>	<b>USER CVS</b>	
	# of Repetitions ACTIVE	<i>0</i>	User CV6	<i>1.00</i>
<b>SAT</b>		<b>MULTI-PHASE</b>		
Tag Type	<i>None</i>	# of Phases	<i>3</i>	
<b>TRICKS</b>		Seperate Series	<i>1</i>	
Pause On/Off	<i>On</i>	Trigger Delay without AV	<i>6</i>	
Auto Subtract	<i>4099</i>	Mask Phase	<i>1</i>	
Auto SCIC	<i>2</i>	Mask Pause	<i>1</i>	
		<b>DIFFUSION</b>		
		Recon All Images	<i>On</i>	
		<b>CONTRAST</b>		
		Contrast Yes/No	<i>No</i>	

3D Ax LAVA Flex Pre/AV

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, ARC, Flex</i>
Plane	<i>CORONAL</i>	IDEAL	<i>3</i>
Series Description	<i>3D Cor LAVA Flex BH</i>	SCAN RANGE	
SCAN TIMING		FOV	<i>44.0</i>
Flip Angle	<i>10</i>	Slice Thickness	<i>4.4</i>
TE	<i>Minimum</i>	Location per Slab	<i>56</i>
Number of Echoes	<i>2</i>	Overlap Locations	<i>0</i>
Receiver Bandwidth	<i>142.86</i>	ACQ TIMING	
IMAGE ENHANCE		Freq	<i>300</i>
Filter Choice	<i>A</i>	Phase	<i>200</i>
GATING/TRIGGER		Freq DIR	<i>S/I</i>
Auto Trigger Type	<i>Off</i>	NEX	<i>1.00</i>
FMRI		Phase FOV	<i>1.00</i>
PSD Trigger	<i>Internal</i>	Auto Shim	<i>Auto</i>
View Order	<i>Bottom/Up</i>	Phase Correction	<i>No</i>
# of Repetitions REST	<i>0</i>	USER CVS	
# of Repetitions ACTIVE	<i>0</i>	User CV4	<i>6.60</i>
SAT		User CV6	<i>1.00</i>
Tag Type	<i>None</i>	MULTI-PHASE	
TRICKS		Seperate Series	<i>0</i>
Pause On/Off	<i>On</i>	Trigger Delay without AV	<i>6.6</i>
Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>
		DIFFUSION	
		Recon All Images	<i>On</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

3D Cor LAVA Flex BH

3D Cor LAVA Flex BH

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

3D Ax LAVA Flex 3MIN	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		3D Ax LAVA Flex 3MIN
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, ARC, Flex</i>	
	Plane	<i>AXIAL</i>	IDEAL	<i>3</i>	
	Series Description	<i>3D Ax LAVA Flex 3MIN</i>	<b>SCAN RANGE</b>		
	<b>SCAN TIMING</b>		FOV	<i>42.0</i>	
	Flip Angle	<i>10</i>	Slice Thickness	<i>4.4</i>	
	TE	<i>Minimum</i>	Location per Slab	<i>56</i>	
	Number of Echoes	<i>2</i>	Overlap Locations	<i>0</i>	
	Receiver Bandwidth	<i>142.86</i>	<b>ACQ TIMING</b>		
	<b>IMAGE ENHANCE</b>		Freq	<i>300</i>	
	Filter Choice	<i>A</i>	Phase	<i>200</i>	
	<b>GATING/TRIGGER</b>		Freq DIR	<i>R/L</i>	
	Auto Trigger Type	<i>Off</i>	NEX	<i>1.00</i>	
	<b>FMRI</b>		Phase FOV	<i>1.10</i>	
	PSD Trigger	<i>Internal</i>	Auto Shim	<i>Auto</i>	
	View Order	<i>Bottom/Up</i>	Phase Correction	<i>No</i>	
	# of Repetitions REST	<i>0</i>	<b>USER CVS</b>		
	# of Repetitions ACTIVE	<i>0</i>	User CV4	<i>6.60</i>	
	<b>SAT</b>		User CV6	<i>1.00</i>	
	Tag Type	<i>None</i>	<b>MULTI-PHASE</b>		
	<b>TRICKS</b>		Seperate Series	<i>0</i>	
	Pause On/Off	<i>On</i>	Trigger Delay without AV	<i>6.6</i>	
	Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>	
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>		
		<b>DIFFUSION</b>			
		Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

3D Ax LAVA Flex 6MIN	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		3D Ax LAVA Flex 6MIN
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, ARC, Flex</i>	
	Plane	<i>AXIAL</i>	IDEAL	<i>3</i>	
	Series Description	<i>3D Ax LAVA Flex 6MIN</i>	<b>SCAN RANGE</b>		
	<b>SCAN TIMING</b>		FOV	<i>42.0</i>	
	Flip Angle	<i>10</i>	Slice Thickness	<i>4.4</i>	
	TE	<i>Minimum</i>	Location per Slab	<i>56</i>	
	Number of Echoes	<i>2</i>	Overlap Locations	<i>0</i>	
	Receiver Bandwidth	<i>142.86</i>	<b>ACQ TIMING</b>		
	<b>IMAGE ENHANCE</b>		Freq	<i>300</i>	
	Filter Choice	<i>A</i>	Phase	<i>200</i>	
	<b>GATING/TRIGGER</b>		Freq DIR	<i>R/L</i>	
	Auto Trigger Type	<i>Off</i>	NEX	<i>1.00</i>	
	<b>FMRI</b>		Phase FOV	<i>1.10</i>	
	PSD Trigger	<i>Internal</i>	Auto Shim	<i>Auto</i>	
	View Order	<i>Bottom/Up</i>	Phase Correction	<i>No</i>	
	# of Repetitions REST	<i>0</i>	<b>USER CVS</b>		
	# of Repetitions ACTIVE	<i>0</i>	User CV4	<i>6.60</i>	
	<b>SAT</b>		User CV6	<i>1.00</i>	
	Tag Type	<i>None</i>	<b>MULTI-PHASE</b>		
	<b>TRICKS</b>		Seperate Series	<i>0</i>	
	Pause On/Off	<i>On</i>	Trigger Delay without AV	<i>6.6</i>	
	Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>	
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>		
		<b>DIFFUSION</b>			
		Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

3D Ax LAVA Flex 10 MIN	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, ARC, Flex</i>
	Plane	<i>AXIAL</i>	IDEAL	<i>3</i>
	Series Description	<i>3D Ax LAVA Flex 10 MIN</i>	<b>SCAN RANGE</b>	
	<b>SCAN TIMING</b>		FOV	<i>42.0</i>
	Flip Angle	<i>12</i>	Slice Thickness	<i>4.4</i>
	TE	<i>Minimum</i>	Location per Slab	<i>56</i>
	Number of Echoes	<i>2</i>	Overlap Locations	<i>0</i>
	Receiver Bandwidth	<i>142.86</i>	<b>ACQ TIMING</b>	
	<b>IMAGE ENHANCE</b>		Freq	<i>300</i>
	Filter Choice	<i>A</i>	Phase	<i>200</i>
	<b>GATING/TRIGGER</b>		Freq DIR	<i>R/L</i>
	Auto Trigger Type	<i>Off</i>	NEX	<i>1.00</i>
	<b>FMRI</b>		Phase FOV	<i>1.10</i>
	PSD Trigger	<i>Internal</i>	Auto Shim	<i>Auto</i>
	View Order	<i>Bottom/Up</i>	Phase Correction	<i>No</i>
	# of Repetitions REST	<i>0</i>	<b>USER CVS</b>	
	# of Repetitions ACTIVE	<i>0</i>	User CV4	<i>6.60</i>
<b>SAT</b>		User CV6	<i>1.00</i>	
Tag Type	<i>None</i>	<b>MULTI-PHASE</b>		
<b>TRICKS</b>		Seperate Series	<i>0</i>	
Pause On/Off	<i>On</i>	Trigger Delay without AV	<i>6.6</i>	
Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>	
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>	
		<b>DIFFUSION</b>		
		Recon All Images	<i>On</i>	
		<b>CONTRAST</b>		
		Contrast Yes/No	<i>No</i>	

3D Ax LAVA Flex 10 MIN



McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

--OPTIONS--	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, EDR, TRF, Fast, SS, FR, ARC</i>
	Plane	<i>3-PLANE</i>	<b>SCAN RANGE</b>	
	Series Description	<i>--OPTIONS--</i>	FOV	<i>48.0</i>
	<b>SCAN TIMING</b>		Slice Thickness	<i>8.0</i>
	TE	<i>80.0</i>	Slice Spacing	<i>5.0</i>
	Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>	
	TR	<i>525.0</i>	Freq	<i>384</i>
	Receiver Bandwidth	<i>83.33</i>	Phase	<i>160</i>
	<b>IMAGE ENHANCE</b>		Freq DIR	<i>Unswap</i>
	Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>
	<b>GATING/TRIGGER</b>		Phase FOV	<i>1.00</i>
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
	<b>FMRI</b>		Phase Correction	<i>No</i>
	PSD Trigger	<i>Internal</i>	<b>USER CVS</b>	
	View Order	<i>Bottom/Up</i>	User CV1	<i>1.00</i>
	# of Repetitions REST	<i>0</i>	<b>MULTI-PHASE</b>	
	# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
	<b>SAT</b>		Mask Phase	<i>0</i>
	Tag Type	<i>None</i>	Mask Pause	<i>0</i>
	<b>TRICKS</b>		<b>DIFFUSION</b>	
	Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>
	Auto Subtract	<i>0</i>	<b>CONTRAST</b>	
Auto SCIC	<i>Off</i>	Contrast Yes/No	<i>No</i>	

--OPTIONS--

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	DISCO
Coil Configuration	Body	Imaging Options	EDR, Fast, ZIP2, MPhVar, ARC, Flex
Plane	AXIAL	IDEAL	3
Series Description	3D Ax DISCO Dyn Mph 3 BH	SCAN RANGE	
SCAN TIMING		FOV	42.0
Flip Angle	12	Slice Thickness	4.4
TE	Minimum	Location per Slab	56
Number of Echoes	2	Overlap Locations	0
Receiver Bandwidth	166.67	ACQ TIMING	
IMAGE ENHANCE		Freq	320
Filter Choice	A	Phase	192
GATING/TRIGGER		Freq DIR	R/L
Auto Trigger Type	Off	Phase FOV	0.80
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
View Order	Bottom/Up	USER CVS	
# of Repetitions REST	0	User CV11	2.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		# of Phases	5
Tag Type	None	Seperate Series	0
TRICKS		Trigger Delay without AV	0
Pause On/Off	On	Mask Phase	1
Auto Subtract	0	Mask Pause	1
Auto SCIC	2	DIFFUSION	
OTHERS		Recon All Images	On
Protocol Notes	<p>Scan Time Dislay:                      1st scan time                      = total time for wash in and delay phases.                      2nd scan time                      = mask scan time/non wash in phases scan time.                      If "Mask" and "Pause After Mask" are selected the scanner will pause after completeing the mask portion and wait for the scan button to be pressed before starting the wash in portions.                      Phases per Locations:                      = wash in and deayed phases                      (Total number of phases)                      Wash-in-Phases:                      = number of phases for the wash in portion.                      - After the pre/mask is complete the scanner will pause.                      - Press the scan button to begin the wash in portions                      - The delayed phases will start after the prescribed delay time runs down.</p>	CONTRAST	
		Contrast Yes/No	No

3D Ax DISCO Dyn Mph 3 BH

3D Ax DISCO Dyn Mph 3 BH

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

Ax Dual Echo BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Ax Dual Echo BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ARC</i>	
	Plane	<i>AXIAL</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Ax Dual Echo BH</i>	FOV	<i>46.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>7.0</i>	
	Flip Angle	<i>55</i>	Slice Spacing	<i>1.0</i>	
	Number of Echoes	<i>2</i>	<b>ACQ TIMING</b>		
	TR	<i>150.0</i>	Freq	<i>288</i>	
	<b>IMAGE ENHANCE</b>		Phase	<i>192</i>	
	Filter Choice	<i>A</i>	Freq DIR	<i>R/L</i>	
	<b>GATING/TRIGGER</b>		Phase FOV	<i>0.80</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>	
	<b>MULTI-PHASE</b>		Phase Correction	<i>No</i>	
	Seperate Series	<i>0</i>	<b>FMRI</b>		
	Mask Phase	<i>0</i>	PSD Trigger	<i>Internal</i>	
	Mask Pause	<i>0</i>	View Order	<i>Bottom/Up</i>	
	<b>DIFFUSION</b>		# of Repetitions REST	<i>0</i>	
	Recon All Images	<i>On</i>	# of Repetitions ACTIVE	<i>0</i>	
<b>CONTRAST</b>		<b>SAT</b>			
Contrast Yes/No	<i>No</i>	Tag Type	<i>None</i>		
		<b>TRICKS</b>			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>2</i>		

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FRFSE-XL</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ZIP2, T2P, Asset, FR, Nav, MRCP</i>
Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
Series Description	<i>3D Obi MRCP Nav</i>	FOV	<i>34.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>1.4</i>
TE	<i>900.0</i>	Location per Slab	<i>64</i>
Number of Echoes	<i>1</i>	Overlap Locations	<i>0</i>
Echo Train Length	<i>120</i>	<b>ACQ TIMING</b>	
Receiver Bandwidth	<i>41.67</i>	Freq	<i>288</i>
<b>IMAGE ENHANCE</b>		Phase	<i>288</i>
Filter Choice	<i>None</i>	Freq DIR	<i>Unswap</i>
<b>GATING/TRIGGER</b>		NEX	<i>1.00</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>
Resp. Trigger Window	<i>30</i>	Auto Shim	<i>Auto</i>
<b>MULTI-PHASE</b>		Phase Correction	<i>No</i>
Seperate Series	<i>0</i>	<b>FMRI</b>	
Mask Phase	<i>0</i>	PSD Trigger	<i>Internal</i>
Mask Pause	<i>0</i>	View Order	<i>Bottom/Up</i>
<b>DIFFUSION</b>		# of Repetitions REST	<i>0</i>
Recon All Images	<i>On</i>	# of Repetitions ACTIVE	<i>0</i>
<b>TRACKER</b>		<b>SAT</b>	
Tracker Length	<i>120.0</i>	Tag Type	<i>None</i>
Tracker Thickness	<i>20.0</i>	Fat/Water Saturation	<i>Fat</i>
		<b>TRICKS</b>	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		<b>CONTRAST</b>	
		Contrast Yes/No	<i>No</i>



3D Obi MRCP Nav

3D Obi MRCP Nav



OTHERS

**Protocol Notes**



*Auto Tracker placement is turned on  
Tracker will be placed automatically  
Manual Placement  
Instructions:  
Navigator Tracker  
Placement Considerations  
\* Place navigator tracker on a breath hold localizer for the optimum performance\*  
- Center of tracker should be 1.5 cm below diaphragm if localizer was acquired with expiration breath hold  
- Center of tracker should be 1.5 cm above diaphragm if localizer was acquired with inspiration breath hold  
- Center of tracker should be placed on the diaphragm if localizer was acquired free breathing  
\*Please review operator manual for additional information  
Expiration-Breath Hold Tracker Placement*





*Inspiration-Breath Hold Tracker Placement*




*Free Breathing Tracker Placement*








## Easy Troubleshooting fo



-> No need to change anything since green waveform reflects respiratory status



-> Increase Track Thickness

McLeod Pavilion

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

Radial SSFSE ASPIR Thick BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Fast, SS, Asset</i>
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
	Series Description	<i>Radial SSFSE ASPIR Thick BH</i>	FOV	<i>38.0</i>
	<b>SCAN TIMING</b>		Slice Thickness	<i>40.0</i>
	TE	<i>900.0</i>	Slice Spacing	<i>0.0</i>
	Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>	
	TR	<i>Minimum</i>	Freq	<i>320</i>
	TI	<i>170</i>	Phase	<i>256</i>
	Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>Unswap</i>
	<b>IMAGE ENHANCE</b>		NEX	<i>1.00</i>
	Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>1</i>
	<b>GATING/TRIGGER</b>		Phase FOV	<i>1.00</i>
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
	<b>FMRI</b>		Phase Correction	<i>No</i>
	PSD Trigger	<i>Internal</i>	<b>USER CVS</b>	
	Slice Order	<i>Interleaved</i>	User CV1	<i>1.00</i>
	View Order	<i>Bottom/Up</i>	User CV2	<i>264.00</i>
# of Repetitions REST	<i>0</i>	User CV13	<i>1.00</i>	
# of Repetitions ACTIVE	<i>0</i>	<b>MULTI-PHASE</b>		
<b>SAT</b>		Seperate Series	<i>0</i>	
Tag Type	<i>None</i>	Mask Phase	<i>0</i>	
Fat/Water Saturation	<i>Fat Special</i>	Mask Pause	<i>0</i>	
<b>TRICKS</b>		<b>DIFFUSION</b>		
Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>	
Auto Subtract	<i>0</i>	<b>CONTRAST</b>		
Auto SCIC	<i>Off</i>	Contrast Yes/No	<i>No</i>	

Radial SSFSE ASPIR Thick BH

Protocol: adult\_abdomen\_LIVER BH 4/3/2018

3D Ax FAT FRACTION IDEAL IQ BH	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Feet First	Imaging Mode	3D
	Patient Position	Supine	Pulse Sequence	IDEAL IQ
	Coil Configuration	Body	Imaging Options	EDR, Fast, IDEAL, ARC
	Plane	AXIAL	IDEAL	195
	Series Description	3D Ax FAT FRACTION IDEAL IQ BH	<b>SCAN RANGE</b>	
	<b>SCAN TIMING</b>		FOV	44.0
	Flip Angle	3	Slice Thickness	10.0
	TE	Min Full	Location per Slab	32
	Number of Echoes	6	Overlap Locations	0
	Echo Train Length	3	<b>ACQ TIMING</b>	
	Number of Shots	2	Freq	160
	Receiver Bandwidth	111.11	Phase	160
	<b>IMAGE ENHANCE</b>		Freq DIR	R/L
	Filter Choice	None	NEX	0.50
	<b>GATING/TRIGGER</b>		Phase FOV	0.90
	Auto Trigger Type	Off	Auto Shim	Auto
	<b>FMRI</b>		Phase Correction	No
	PSD Trigger	Internal	<b>USER CVS</b>	
	Slice Order	Interleaved	User CV4	5.60
View Order	Bottom/Up	<b>MULTI-PHASE</b>		
# of Repetitions REST	0	Seperate Series	0	
# of Repetitions ACTIVE	0	Trigger Delay without AV	5.6	
<b>SAT</b>		Mask Phase	0	
Tag Type	None	Mask Pause	0	
<b>TRICKS</b>		<b>DIFFUSION</b>		
Pause On/Off	On	Recon All Images	On	
Auto Subtract	0	<b>CONTRAST</b>		
Auto SCIC	Off	Contrast Yes/No	No	
<b>OTHERS</b>				
Protocol Notes	<p>* Number of echoes = Number of shots x ETL The "Auto" Flip angle feature optimizes this sequence for maximal SNR, minimal T1 weighting and maximum fat fraction accuracy for the liver. If this sequence is used outside of the liver, the fixed flip angle value of 3 degrees is recommended otherwise fat fraction values may be overestimated R2* value may be overestimated after injection of contrast media. *Fat fraction value may change after injection of contrast media especially when the flip angle is sufficiently high. Please use auto flip angle option</p>			

3D Ax FAT FRACTION IDEAL IQ BH