

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

3Plane Loc FIESTA	PATIENT POSITION		IMAGING PARAMETERS		3Plane Loc FIESTA
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, Fast</i>	
	Plane	<i>3-PLANE</i>	SCAN RANGE		
	Series Description	<i>3Plane Loc FIESTA</i>	FOV	<i>48.0</i>	
	SCAN TIMING		Slice Thickness	<i>8.0</i>	
	Flip Angle	<i>45</i>	Slice Spacing	<i>5.0</i>	
	TE	<i>Minimum</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>160</i>	
	Receiver Bandwidth	<i>100.00</i>	Phase	<i>140</i>	
	IMAGE ENHANCE		Freq DIR	<i>Unswap</i>	
	Filter Choice	<i>None</i>	NEX	<i>1.00</i>	
	GATING/TRIGGER		# of Acq. Before Pause	<i>0</i>	
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>	
	MULTI-PHASE		Auto Shim	<i>Auto</i>	
	Seperate Series	<i>0</i>	Phase Correction	<i>No</i>	
	Mask Phase	<i>0</i>	FMRI		
	Mask Pause	<i>0</i>	PSD Trigger	<i>Internal</i>	
	DIFFUSION		View Order	<i>Bottom/Up</i>	
Recon All Images	<i>On</i>	# of Repetitions REST	<i>0</i>		
CONTRAST		# of Repetitions ACTIVE	<i>0</i>		
Contrast Yes/No	<i>No</i>	SAT			
		Tag Type	<i>None</i>		
		TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, EDR, Fast, ZIP512, Asset</i>
Plane	<i>AXIAL</i>	SCAN RANGE	
Series Description	<i>Ax FIESTA BH Loc</i>	FOV	<i>36.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>55</i>	Slice Spacing	<i>0.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>132</i>
Receiver Bandwidth	<i>100.00</i>	Phase	<i>132</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>0</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>
FMRI		Auto Shim	<i>Auto</i>
PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>
View Order	<i>Bottom/Up</i>	USER CVS	
# of Repetitions REST	<i>0</i>	User CV20	<i>1.00</i>
# of Repetitions ACTIVE	<i>0</i>	MULTI-PHASE	
SAT		Seperate Series	<i>0</i>
Tag Type	<i>None</i>	Mask Phase	<i>0</i>
TRICKS		Mask Pause	<i>0</i>
Pause On/Off	<i>On</i>	DIFFUSION	
Auto Subtract	<i>0</i>	Optimized TE	<i>Yes</i>
Auto SCIC	<i>2</i>	Recon All Images	<i>On</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

AX FIESTA BH Loc

AX FIESTA BH Loc

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

2Ch LOC BH	PATIENT POSITION		IMAGING PARAMETERS		2Ch LOC BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>2Ch LOC BH</i>	FOV	<i>41.0</i>	
	SCAN TIMING		Slice Thickness	<i>8.0</i>	
	Flip Angle	<i>50</i>	Slice Spacing	<i>0.0</i>	
	TE	<i>Minimum</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>192</i>	
	Receiver Bandwidth	<i>100.00</i>	Phase	<i>224</i>	
	IMAGE ENHANCE		Freq DIR	<i>S/I</i>	
	Filter Choice	<i>None</i>	NEX	<i>1.00</i>	
	GATING/TRIGGER		Phase FOV	<i>1.00</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>On</i>	
	Trigger Type	<i>1</i>	Phase Correction	<i>No</i>	
	Views per Segment	<i>20</i>	FMRI		
	MULTI-PHASE		PSD Trigger	<i>Internal</i>	
	Seperate Series	<i>0</i>	Slice Order	<i>Interleaved</i>	
	Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>	
Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>		
DIFFUSION		# of Repetitions ACTIVE	<i>0</i>		
Recon All Images	<i>On</i>	SAT			
CINE		Tag Type	<i>None</i>		
# of Cardiac Phases to Reconstruct	<i>20</i>	TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		
		CONTRAST			
		Contrast Yes/No	<i>No</i>		

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SEMI 4CH LOC BH</i>	FOV	<i>41.0</i>
SCAN TIMING		Slice Thickness	<i>7.0</i>
Flip Angle	<i>50</i>	Slice Spacing	<i>0.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>192</i>
Receiver Bandwidth	<i>100.00</i>	Phase	<i>224</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		Phase FOV	<i>0.90</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>On</i>
Trigger Type	<i>1</i>	Phase Correction	<i>No</i>
Views per Segment	<i>20</i>	FMRI	
MULTI-PHASE		PSD Trigger	<i>Internal</i>
Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
DIFFUSION		SAT	
Recon All Images	<i>On</i>	Tag Type	<i>None</i>
CINE		TRICKS	
# of Cardiac Phases to Reconstruct	<i>20</i>	Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

SEMI 4CH LOC BH

SEMI 4CH LOC BH

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SA LOC BH</i>	FOV	<i>41.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>50</i>	Slice Spacing	<i>0.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>160</i>
Receiver Bandwidth	<i>100.00</i>	Phase	<i>160</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		Phase FOV	<i>0.80</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>On</i>
Trigger Type	<i>1</i>	Phase Correction	<i>No</i>
Views per Segment	<i>20</i>	FMRI	
MULTI-PHASE		PSD Trigger	<i>Internal</i>
Seperate Series	<i>0</i>	Slice Order	<i>Interleaved</i>
Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>
DIFFUSION		# of Repetitions ACTIVE	<i>0</i>
Recon All Images	<i>On</i>	SAT	
CINE		Tag Type	<i>None</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

SA LOC BH

SA LOC BH

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

4CH LOC BH	PATIENT POSITION		IMAGING PARAMETERS		4CH LOC BH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>4CH LOC BH</i>	FOV	<i>41.0</i>	
	SCAN TIMING		Slice Thickness	<i>7.0</i>	
	Flip Angle	<i>50</i>	Slice Spacing	<i>0.0</i>	
	TE	<i>Minimum</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>192</i>	
	Receiver Bandwidth	<i>100.00</i>	Phase	<i>224</i>	
	IMAGE ENHANCE		Freq DIR	<i>R/L</i>	
	Filter Choice	<i>None</i>	NEX	<i>1.00</i>	
	GATING/TRIGGER		Phase FOV	<i>0.90</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>On</i>	
	Trigger Type	<i>1</i>	Phase Correction	<i>No</i>	
	Views per Segment	<i>20</i>	FMRI		
	MULTI-PHASE		PSD Trigger	<i>Internal</i>	
	Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>	
	Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>	
	Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>	
	DIFFUSION		SAT		
	Recon All Images	<i>On</i>	Tag Type	<i>None</i>	
	CINE		TRICKS		
	# of Cardiac Phases to Reconstruct	<i>20</i>	Pause On/Off	<i>On</i>	
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		
		CONTRAST			
		Contrast Yes/No	<i>No</i>		

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

SA LOC FREE BREATH	PATIENT POSITION		IMAGING PARAMETERS		SA LOC FREE BREATH
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, Fast, Asset</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>SA LOC FREE BREATH</i>	FOV	<i>41.0</i>	
	SCAN TIMING		Slice Thickness	<i>8.0</i>	
	Flip Angle	<i>50</i>	Slice Spacing	<i>0.0</i>	
	TE	<i>Minimum</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>160</i>	
	Receiver Bandwidth	<i>100.00</i>	Phase	<i>160</i>	
	IMAGE ENHANCE		Freq DIR	<i>S/I</i>	
	Filter Choice	<i>None</i>	NEX	<i>1.00</i>	
	GATING/TRIGGER		# of Acq. Before Pause	<i>0</i>	
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.80</i>	
	MULTI-PHASE		Auto Shim	<i>On</i>	
	Seperate Series	<i>0</i>	Phase Correction	<i>No</i>	
	Mask Phase	<i>0</i>	FMRI		
	Mask Pause	<i>0</i>	PSD Trigger	<i>Internal</i>	
	DIFFUSION		View Order	<i>Bottom/Up</i>	
Recon All Images	<i>On</i>	# of Repetitions REST	<i>0</i>		
CONTRAST		# of Repetitions ACTIVE	<i>0</i>		
Contrast Yes/No	<i>No</i>	SAT			
		Tag Type	<i>None</i>		
		TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SA FIESTA CINE BH</i>	FOV	<i>44.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>45</i>	Slice Spacing	<i>1.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>192</i>
Receiver Bandwidth	<i>125.00</i>	Phase	<i>256</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>1</i>
Auto Trigger Type	<i>On</i>	Phase FOV	<i>0.80</i>
Trigger Type	<i>0</i>	Auto Shim	<i>On</i>
Views per Segment	<i>20</i>	Phase Correction	<i>No</i>
MULTI-PHASE		FMRI	
Seperate Series	<i>0</i>	PSD Trigger	<i>Internal</i>
Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>
DIFFUSION		# of Repetitions ACTIVE	<i>0</i>
Recon All Images	<i>On</i>	SAT	
CINE		Tag Type	<i>None</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

SA FIESTA CINE BH

SA FIESTA CINE BH

McLeod Pavilion

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>4CH FIESTA BH</i>	FOV	<i>42.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>45</i>	Slice Spacing	<i>0.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>192</i>
Receiver Bandwidth	<i>142.86</i>	Phase	<i>192</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>1</i>
Auto Trigger Type	<i>On</i>	Phase FOV	<i>0.90</i>
Trigger Type	<i>0</i>	Auto Shim	<i>On</i>
Views per Segment	<i>20</i>	Phase Correction	<i>No</i>
MULTI-PHASE		FMRI	
Seperate Series	<i>0</i>	PSD Trigger	<i>Internal</i>
Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>
DIFFUSION		# of Repetitions ACTIVE	<i>0</i>
Recon All Images	<i>On</i>	SAT	
CINE		Tag Type	<i>None</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

4CH FIESTA BH

4CH FIESTA BH

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, Asset</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>3CH FIESTA BH</i>	FOV	<i>38.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>45</i>	Slice Spacing	<i>0.0</i>
TE	<i>Minimum</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>192</i>
Receiver Bandwidth	<i>142.86</i>	Phase	<i>192</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>1</i>
Auto Trigger Type	<i>On</i>	Phase FOV	<i>1.00</i>
Trigger Type	<i>0</i>	Auto Shim	<i>On</i>
Views per Segment	<i>20</i>	Phase Correction	<i>No</i>
MULTI-PHASE		FMRI	
Seperate Series	<i>0</i>	PSD Trigger	<i>Internal</i>
Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>
DIFFUSION		# of Repetitions ACTIVE	<i>0</i>
Recon All Images	<i>On</i>	SAT	
CINE		Tag Type	<i>None</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

3CH FIESTA BH

3CH FIESTA BH

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, ZIP512, SS, BSP, Asset</i>
Plane	<i>AXIAL</i>	SCAN RANGE	
Series Description	<i>AX T2 FB BLACK BLOOD SSFSE</i>	FOV	<i>36.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
TE	<i>50.0</i>	Slice Spacing	<i>0.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
Blood Suppression TI	<i>737</i>	Freq	<i>212</i>
Receiver Bandwidth	<i>83.33</i>	Phase	<i>150</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>A</i>	# of Acq. Before Pause	<i>0</i>
GATING/TRIGGER		Phase FOV	<i>0.90</i>
Auto Trigger Type	<i>On</i>	Auto Shim	<i>Auto</i>
Trigger Type	<i>0</i>	Phase Correction	<i>No</i>
Minimum Trigger Delay	<i>2</i>	USER CVS	
Cardiac Phase	<i>1</i>	User CV2	<i>240.00</i>
Cardiac Slices	<i>1</i>	User CV13	<i>1.00</i>
FMRI		MULTI-PHASE	
PSD Trigger	<i>Internal</i>	Seperate Series	<i>0</i>
View Order	<i>Bottom/Up</i>	Mask Phase	<i>0</i>
# of Repetitions REST	<i>0</i>	Mask Pause	<i>0</i>
# of Repetitions ACTIVE	<i>0</i>	DIFFUSION	
SAT		Recon All Images	<i>On</i>
Tag Type	<i>None</i>	CONTRAST	
TRICKS		Contrast Yes/No	<i>No</i>
# of RR Intervals	<i>2</i>		
Trigger Delay	<i>748</i>		
Pause On/Off	<i>On</i>		
Auto Subtract	<i>0</i>		
Auto SCIC	<i>2</i>		

AX T2 FB BLACK BLOOD SSFSE

AX T2 FB BLACK BLOOD SSFSE

Protocol: adult_chest_CARDIAC MASS

AX Triple IR FB BLACK BLOOD	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, ZIP512, BSP, Asset, IrP</i>
	Plane	<i>AXIAL</i>	PSD Name	<i>FSE-XL</i>
	Series Description	<i>AX Triple IR FB BLACK BLOOD</i>	SCAN RANGE	
	SCAN TIMING		FOV	<i>36.0</i>
	Flip Angle	<i>107</i>	Slice Thickness	<i>8.0</i>
	TE	<i>65.0</i>	Slice Spacing	<i>0.0</i>
	Number of Echoes	<i>1</i>	ACQ TIMING	
	TI	<i>220</i>	Freq	<i>224</i>
	Echo Train Length	<i>23</i>	Phase	<i>192</i>
	Blood Suppression TI	<i>1053</i>	Freq DIR	<i>R/L</i>
	Receiver Bandwidth	<i>50.00</i>	Fat Shift DIR	<i>Normal (L)</i>
	IMAGE ENHANCE		NEX	<i>1.00</i>
	Filter Choice	<i>A</i>	# of Acq. Before Pause	<i>0</i>
	GATING/TRIGGER		Phase FOV	<i>0.90</i>
	Auto Trigger Type	<i>On</i>	Auto Shim	<i>Auto</i>
	Trigger Type	<i>0</i>	Phase Correction	<i>No</i>
	Minimum Trigger Delay	<i>2</i>	USER CVS	
	Cardiac Phase	<i>1</i>	User CV14	<i>3.33</i>
	Cardiac Slices	<i>1</i>	User CV22	<i>1.00</i>
	FMRI		MULTI-PHASE	
	PSD Trigger	<i>Internal</i>	Seperate Series	<i>0</i>
	View Order	<i>Bottom/Up</i>	Mask Phase	<i>0</i>
# of Repetitions REST	<i>0</i>	Mask Pause	<i>0</i>	
# of Repetitions ACTIVE	<i>0</i>	DIFFUSION		
SAT		Recon All Images	<i>On</i>	
Tag Type	<i>None</i>	CONTRAST		
TRICKS		Contrast Yes/No	<i>No</i>	
# of RR Intervals	<i>2</i>			
Trigger Delay	<i>1065</i>			
Pause On/Off	<i>On</i>			
Auto Subtract	<i>0</i>			
Auto SCIC	<i>2</i>			

AX Triple IR FB BLACK BLOOD

Protocol: adult_chest_CARDIAC MASS

4CH SS PS MDE PRE	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
	Plane	<i>OBLIQUE</i>	SCAN RANGE	
	Series Description	<i>4CH SS PS MDE PRE</i>	FOV	<i>36.0</i>
	SCAN TIMING		Slice Thickness	<i>8.0</i>
	Flip Angle	<i>20</i>	Slice Spacing	<i>0.0</i>
	Number of Echoes	<i>1</i>	ACQ TIMING	
	TI	<i>600</i>	Freq	<i>156</i>
	Receiver Bandwidth	<i>62.50</i>	Phase	<i>140</i>
	IMAGE ENHANCE		Freq DIR	<i>R/L</i>
	Filter Choice	<i>A</i>	NEX	<i>1.00</i>
	GATING/TRIGGER		# of Acq. Before Pause	<i>4</i>
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.80</i>
	Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
	Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
	Cardiac Phase	<i>1</i>	FMRI	
	Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
	MULTI-PHASE		View Order	<i>Bottom/Up</i>
	Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
	Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
	Mask Pause	<i>0</i>	SAT	
	DIFFUSION		Tag Type	<i>None</i>
	Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>	
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>605</i>	
		Pause On/Off	<i>On</i>	
		Auto Subtract	<i>0</i>	
		Auto SCIC	<i>Off</i>	
		CONTRAST		
		Contrast Yes/No	<i>No</i>	

4CH SS PS MDE PRE

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, MPh, Asset, SrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SA PERFUSION FB SCOUT</i>	FOV	<i>36.0</i>
SCAN TIMING		Slice Thickness	<i>10.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>1.0</i>
TE	<i>Min Full</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>144</i>
TI	<i>150</i>	Phase	<i>132</i>
Receiver Bandwidth	<i>62.50</i>	Freq DIR	<i>S/I</i>
IMAGE ENHANCE		NEX	<i>0.75</i>
Filter Choice	<i>None</i>	Phase FOV	<i>0.80</i>
GATING/TRIGGER		Auto Shim	<i>Auto</i>
Auto Trigger Type	<i>Off</i>	Phase Correction	<i>No</i>
Trigger Type	<i>1</i>	USER CVS	
Minimum Trigger Delay	<i>2</i>	User CV16	<i>1.00</i>
Cardiac Phase	<i>1</i>	MULTI-PHASE	
Cardiac Slices	<i>1</i>	Slice per Location	<i>1</i>
FMRI		Seperate Series	<i>0</i>
PSD Trigger	<i>Internal</i>	Mask Phase	<i>0</i>
Slice Order	<i>Interleaved</i>	Mask Pause	<i>0</i>
View Order	<i>Bottom/Up</i>	DIFFUSION	
# of Repetitions REST	<i>0</i>	Recon All Images	<i>On</i>
# of Repetitions ACTIVE	<i>0</i>	CONTRAST	
SAT		Contrast Yes/No	<i>No</i>
Tag Type	<i>None</i>		
TRICKS			
# of RR Intervals	<i>2</i>		
Trigger Delay	<i>10</i>		
Pause On/Off	<i>On</i>		
Auto Subtract	<i>0</i>		
Auto SCIC	<i>Off</i>		

SA PERFUSION FB SCOUT

SA PERFUSION FB SCOUT

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, MPh, Asset, SrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SA PERFUSION BH +C</i>	FOV	<i>36.0</i>
SCAN TIMING		Slice Thickness	<i>10.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>4.0</i>
TE	<i>Min Full</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>144</i>
TI	<i>150</i>	Phase	<i>132</i>
Receiver Bandwidth	<i>62.50</i>	Freq DIR	<i>S/I</i>
IMAGE ENHANCE		NEX	<i>0.75</i>
Filter Choice	<i>None</i>	Phase FOV	<i>0.80</i>
GATING/TRIGGER		Auto Shim	<i>Auto</i>
Auto Trigger Type	<i>Off</i>	Phase Correction	<i>No</i>
Trigger Type	<i>1</i>	USER CVS	
Minimum Trigger Delay	<i>2</i>	User CV16	<i>1.00</i>
Cardiac Phase	<i>1</i>	MULTI-PHASE	
Cardiac Slices	<i>1</i>	Slice per Location	<i>30</i>
FMRI		Seperate Series	<i>0</i>
PSD Trigger	<i>Internal</i>	Mask Phase	<i>0</i>
View Order	<i>Bottom/Up</i>	Mask Pause	<i>0</i>
# of Repetitions REST	<i>0</i>	DIFFUSION	
# of Repetitions ACTIVE	<i>0</i>	Recon All Images	<i>On</i>
SAT		CONTRAST	
Tag Type	<i>None</i>	Contrast Yes/No	<i>Yes</i>
TRICKS			
# of RR Intervals	<i>2</i>		
Trigger Delay	<i>10</i>		
Pause On/Off	<i>On</i>		
Auto Subtract	<i>0</i>		
Auto SCIC	<i>Off</i>		

SA PERFUSION BH +C

SA PERFUSION BH +C

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>4CH SS PS MDE BH+C</i>	FOV	<i>36.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>0.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TI	<i>600</i>	Freq	<i>156</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>140</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>4</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.80</i>
Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
MULTI-PHASE		View Order	<i>Bottom/Up</i>
Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Pause	<i>0</i>	SAT	
DIFFUSION		Tag Type	<i>None</i>
Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>605</i>
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

4CH SS PS MDE BH+C

4CH SS PS MDE BH+C

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>LVOT SS PS MDE BH +C</i>	FOV	<i>42.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>0.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TI	<i>600</i>	Freq	<i>156</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>140</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>0</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.90</i>
Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
MULTI-PHASE		View Order	<i>Bottom/Up</i>
Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Pause	<i>0</i>	SAT	
DIFFUSION		Tag Type	<i>None</i>
Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>605</i>
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

LVOT SS PS MDE BH +C

LVOT SS PS MDE BH +C

Protocol: adult_chest_CARDIAC MASS

Ax 4DFLOW postCE FB	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>4DFLOW</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, ZIP2, ARC</i>
	Plane	<i>AXIAL</i>	SCAN RANGE	
	Series Description	<i>Ax 4DFLOW postCE FB</i>	FOV	<i>46.0</i>
	SCAN TIMING		Slice Thickness	<i>2.4</i>
	Flip Angle	<i>15</i>	Location per Slab	<i>88</i>
	TE	<i>Min Full</i>	Overlap Locations	<i>0</i>
	Number of Echoes	<i>1</i>	ACQ TIMING	
	Receiver Bandwidth	<i>62.50</i>	Freq	<i>256</i>
	IMAGE ENHANCE		Phase	<i>256</i>
	Filter Choice	<i>None</i>	Freq DIR	<i>R/L</i>
	GATING/TRIGGER		NEX	<i>4.00</i>
	Auto Trigger Type	<i>On</i>	Phase FOV	<i>0.65</i>
	Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
	Views per Segment	<i>3</i>	Phase Correction	<i>No</i>
	# of Acq Phases	<i>16</i>	USER CVS	
	FMRI		User CV3	<i>1.00</i>
	PSD Trigger	<i>Internal</i>	User CV5	<i>65.00</i>
	View Order	<i>Bottom/Up</i>	User CV6	<i>2.00</i>
	# of Repetitions REST	<i>0</i>	User CV16	<i>10.00</i>
	# of Repetitions ACTIVE	<i>0</i>	MULTI-PHASE	
	SAT		Seperate Series	<i>0</i>
	Tag Type	<i>None</i>	Mask Phase	<i>0</i>
	DIFFUSION		Mask Pause	<i>0</i>
	Recon All Images	<i>On</i>	VASCULAR	
	TRICKS		Flow Analysis	<i>All</i>
Pause On/Off	<i>On</i>	Velocity Encoding	<i>300.0</i>	
Auto Subtract	<i>0</i>	ASSET		
Auto SCIC	<i>Off</i>	Hyperkat Value	<i>8</i>	
CONTRAST		CINE		
Contrast Yes/No	<i>Yes</i>	# of Cardiac Phases to Reconstruct	<i>20</i>	

Ax 4DFLOW postCE FB

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>CINE IR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, Seq, EDR, Fast, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>CINE IR 2RR BH</i>	FOV	<i>40.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>5</i>	Slice Spacing	<i>0.0</i>
TE	<i>Min Full</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>128</i>
Receiver Bandwidth	<i>50.00</i>	Phase	<i>128</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>None</i>	NEX	<i>1.00</i>
GATING/TRIGGER		Phase FOV	<i>0.80</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
Trigger Type	<i>1</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
Views per Segment	<i>14</i>	View Order	<i>Bottom/Up</i>
MULTI-PHASE		# of Repetitions REST	<i>0</i>
Seperate Series	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Phase	<i>0</i>	SAT	
Mask Pause	<i>0</i>	Tag Type	<i>None</i>
DIFFUSION		TRICKS	
Recon All Images	<i>On</i>	# of RR Intervals	<i>2</i>
CINE		Pause On/Off	<i>On</i>
# of Cardiac Phases to Reconstruct	<i>30</i>	Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

CINE IR 2RR BH

CINE IR 2RR BH

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>SA SS PS MDE BH</i>	FOV	<i>44.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>1.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TI	<i>400</i>	Freq	<i>200</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>200</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>4</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.80</i>
Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
MULTI-PHASE		View Order	<i>Bottom/Up</i>
Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Pause	<i>0</i>	SAT	
DIFFUSION		Tag Type	<i>None</i>
Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>405</i>
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

SA SS PS MDE BH

SA SS PS MDE BH

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>2 CH SS PS MDE SPGR</i>	FOV	<i>41.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>0.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TI	<i>400</i>	Freq	<i>156</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>140</i>
IMAGE ENHANCE		Freq DIR	<i>S/I</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>6</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>
Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
MULTI-PHASE		View Order	<i>Bottom/Up</i>
Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Pause	<i>0</i>	SAT	
DIFFUSION		Tag Type	<i>None</i>
Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>405</i>
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

2 CH SS PS MDE SPGR

2 CH SS PS MDE SPGR

Protocol: adult_chest_CARDIAC MASS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>Gat, EDR, Fast, PS, Asset, IrP</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>4CH SS PS MDE SPGR</i>	FOV	<i>38.0</i>
SCAN TIMING		Slice Thickness	<i>8.0</i>
Flip Angle	<i>20</i>	Slice Spacing	<i>0.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TI	<i>425</i>	Freq	<i>156</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>140</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>6</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.90</i>
Trigger Type	<i>1</i>	Auto Shim	<i>Auto</i>
Minimum Trigger Delay	<i>2</i>	Phase Correction	<i>No</i>
Cardiac Phase	<i>1</i>	FMRI	
Cardiac Slices	<i>1</i>	PSD Trigger	<i>Internal</i>
MULTI-PHASE		View Order	<i>Bottom/Up</i>
Seperate Series	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Phase	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
Mask Pause	<i>0</i>	SAT	
DIFFUSION		Tag Type	<i>None</i>
Recon All Images	<i>On</i>	TRICKS	
CINE		# of RR Intervals	<i>2</i>
# of Cardiac Phases to Reconstruct	<i>20</i>	Trigger Delay	<i>430</i>
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>
		CONTRAST	
		Contrast Yes/No	<i>Yes</i>

4CH SS PS MDE SPGR

4CH SS PS MDE SPGR