

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	48Ch Head Coil	Imaging Options	Seq, EDR, TRF, Fast, SS, ARC
Plane	3-PLANE	SCAN RANGE	
Series Description	3-Plane Localizer	FOV	30.0
SCAN TIMING		Slice Thickness	10.0
TE	80.0	Slice Spacing	5.0
Number of Echoes	1	ACQ TIMING	
TR	1200.0	Freq	256
Receiver Bandwidth	83.33	Phase	128
IMAGE ENHANCE		Freq DIR	Unswap
Filter Choice	None	# of Acq. Before Pause	0
GATING/TRIGGER		Phase FOV	1.00
Auto Trigger Type	Off	Auto Shim	Auto
FMRI		Phase Correction	No
PSD Trigger	Internal	USER CVS	
View Order	Bottom/Up	User CV1	1.00
# of Repetitions REST	0	MULTI-PHASE	
# of Repetitions ACTIVE	0	Seperate Series	0
SAT		Mask Phase	0
Tag Type	None	Mask Pause	0
TRICKS		DIFFUSION	
Pause On/Off	On	Recon All Images	On
Auto Subtract	0	CONTRAST	
Auto SCIC	Off	Contrast Yes/No	No

3-Plane Localizer

3-Plane Localizer

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, EPI, DIFF, Asset
Plane	OBLIQUE	SCAN RANGE	
Series Description	Ax DWI ALL B-1000	FOV	24.0
SCAN TIMING		Slice Thickness	4.0
TE	Minimum	Slice Spacing	0.4
Number of Echoes	1	ACQ TIMING	
TR	7424.0	Freq	120
Number of Shots	1	Phase	160
IMAGE ENHANCE		Freq DIR	Swap
Filter Choice	None	Phase FOV	1.00
GATING/TRIGGER		Auto Shim	Auto
Auto Trigger Type	Off	Phase Correction	Yes
FMRI		USER CVS	
PSD Trigger	Internal	User CV0	1.00
View Order	Bottom/Up	User CV5	1.00
# of Repetitions REST	0	User CV17	1.00
# of Repetitions ACTIVE	0	TR Min	2000.0
SAT		TR Max	17000.0
Tag Type	None	MULTI-PHASE	
TRICKS		Seperate Series	0
Pause On/Off	On	Mask Phase	0
Auto Subtract	0	Mask Pause	0
Auto SCIC	2	DIFFUSION	
		Optimized TE	Yes
		Diffusion Directions	All
		Number of Diffusion Directions	3
		Number of T2 Images	1
		Dual Spin Echo	Off
		Recon All Images	On
		CONTRAST	
		Contrast Yes/No	No

Ax DWI ALL B-1000

Ax DWI ALL B-1000

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Gradient Echo
Coil Configuration	48Ch Head Coil	Imaging Options	EDR
Plane	OBLIQUE	SCAN RANGE	
Series Description	Ax T2* GRE	FOV	24.0
SCAN TIMING		Slice Thickness	4.0
Flip Angle	15	Slice Spacing	1.5
TE	Min Full	ACQ TIMING	
Number of Echoes	1	Freq	256
TR	509.0	Phase	160
Receiver Bandwidth	20.83	Freq DIR	A/P
IMAGE ENHANCE		NEX	2.00
Filter Choice	None	Phase FOV	1.00
GATING/TRIGGER		Auto Shim	Auto
Auto Trigger Type	Off	Phase Correction	No
MULTI-PHASE		FMRI	
Seperate Series	0	PSD Trigger	Internal
Mask Phase	0	Slice Order	Interleaved
Mask Pause	0	View Order	Bottom/Up
DIFFUSION		# of Repetitions REST	0
Recon All Images	On	# of Repetitions ACTIVE	0
CONTRAST		SAT	
Contrast Yes/No	No	Tag Type	None
		TRICKS	
		Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	On

Ax T2* GRE

Ax T2* GRE

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	Head;Neck	Imaging Options	FC, EDR, TRF, Fast
Plane	OBLIQUE	SCAN RANGE	
Series Description	Ax T2 FSE	FOV	21.0
SCAN TIMING		Slice Thickness	4.0
Flip Angle	142	Slice Spacing	0.4
TE	102.0	ACQ TIMING	
Number of Echoes	1	Freq	320
TR	5184.0	Phase	320
Echo Train Length	19	Freq DIR	A/P
Receiver Bandwidth	41.67	Fat Shift DIR	Normal (A)
IMAGE ENHANCE		NEX	1.00
Filter Choice	G	Phase FOV	1.00
GATING/TRIGGER		Auto Shim	Auto
Auto Trigger Type	Off	Phase Correction	No
FMRI		Flow Direction Compensation	Slice
PSD Trigger	Internal	USER CVS	
Slice Order	Interleaved	User CV22	1.00
View Order	Bottom/Up	TR Min	3000.0
# of Repetitions REST	0	TR Max	10000.0
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
SAT Location	1	Mask Phase	0
Tag Type	None	Mask Pause	0
TRICKS		DIFFUSION	
Pause On/Off	On	Recon All Images	On
Auto Subtract	0	CONTRAST	
Auto SCIC	2	Contrast Yes/No	No

Ax T2 FSE

Ax T2 FSE

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/w/o

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, TRF, Fast, ARC, T2flair
Plane	OBLIQUE		
Series Description	Ax T2 FLAIR		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	24.0
TE	140.0	Slice Thickness	4.0
Number of Echoes	1	Slice Spacing	1.5
TR	12000.0		
TI	2711	ACQ TIMING	
Echo Train Length	21	Freq	260
Receiver Bandwidth	41.67	Phase	224
IMAGE ENHANCE		Freq DIR	A/P
Filter Choice	A	Fat Shift DIR	Normal (A)
GATING/TRIGGER		NEX	1.00
Auto Trigger Type	Off	Phase FOV	1.00
FMRI		Auto Shim	On
PSD Trigger	Internal	Phase Correction	No
View Order	Bottom/Up	USER CVS	
# of Repetitions REST	0	User CV22	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
		Contrast Yes/No	No

Ax T2 FLAIR

Ax T2 FLAIR

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, TRF, Fast, ARC, T2flair
Plane	OBLIQUE		
Series Description	Cor T2 FLAIR		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	24.0
TE	140.0	Slice Thickness	4.0
Number of Echoes	1	Slice Spacing	1.5
TR	12000.0		
TI	2711	ACQ TIMING	
Echo Train Length	21	Freq	260
Receiver Bandwidth	41.67	Phase	224
IMAGE ENHANCE		Freq DIR	A/P
Filter Choice	A	Fat Shift DIR	Normal (A)
GATING/TRIGGER		NEX	1.00
Auto Trigger Type	Off	Phase FOV	1.00
FMRI		Auto Shim	On
PSD Trigger	Internal	Phase Correction	No
View Order	Bottom/Up	USER CVS	
# of Repetitions REST	0	User CV22	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
		Contrast Yes/No	No

Cor T2 FLAIR

Cor T2 FLAIR

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	BRAVO
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, Fast, ZIP512, ZIP2, ARC, IrP
Plane	AXIAL		
Series Description	3D Ax T1 BRAVO		
SCAN TIMING		SCAN RANGE	
Flip Angle	12	FOV	24.0
Number of Echoes	1	Slice Thickness	1.2
TI	450	Location per Slab	130
Receiver Bandwidth	31.25	Overlap Locations	0
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	G	Freq	256
		Phase	256
GATING/TRIGGER		Freq DIR	A/P
Auto Trigger Type	Off	NEX	1.00
MULTI-PHASE		Phase FOV	1.00
Seperate Series	0	Auto Shim	Auto
Trigger Delay without AV	0	Phase Correction	No
Mask Phase	0	FMRI	
Mask Pause	0	PSD Trigger	Internal
DIFFUSION		View Order	Bottom/Up
Recon All Images	On	# of Repetitions REST	0
CONTRAST		# of Repetitions ACTIVE	0
Contrast Yes/No	No	SAT	
		Tag Type	None
		TRICKS	
		Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	2

3D Ax T1 BRAVO

3D Ax T1 BRAVO

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	48Ch Head Coil	Imaging Options	Seq, EDR, TRF, Fast, SS, ARC
Plane	3-PLANE		
Series Description	--Contrast--		
SCAN TIMING		SCAN RANGE	
TE	80.0	FOV	30.0
Number of Echoes	1	Slice Thickness	10.0
TR	1200.0	Slice Spacing	5.0
Receiver Bandwidth	83.33		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	256
		Phase	128
GATING/TRIGGER		Freq DIR	Unswap
Auto Trigger Type	Off	# of Acq. Before Pause	0
		Phase FOV	1.00
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
View Order	Bottom/Up	USER CVS	
# of Repetitions REST	0	User CV1	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
		Mask Pause	0
TRICKS		DIFFUSION	
Pause On/Off	On	Recon All Images	On
Auto Subtract	0	CONTRAST	
Auto SCIC	Off	Contrast Yes/No	No

McLeod Pavilion

Protocol: adult_head_Seizure Brain w/wo

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	SPGR
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, Fast, ZIP512, ZIP2, Asset, IrP
Plane	OBLIQUE		
Series Description	+ C 3D Ax T1 FSPGR		
SCAN TIMING		SCAN RANGE	
Flip Angle	12	FOV	24.0
TE	Min Full	Slice Thickness	1.2
Number of Echoes	1	Location per Slab	130
TI	450	Overlap Locations	0
Receiver Bandwidth	31.25		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	G	Freq	256
		Phase	256
GATING/TRIGGER		Freq DIR	Unswap
Auto Trigger Type	Off	NEX	1.00
FMRI		Phase FOV	1.00
Initial State	Control	Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
View Order	Bottom/Up	USER CVS	
# of Repetitions REST	0	User CV6	2.00
# of Repetitions ACTIVE	0	User CV23	90.00
DIFFUSION		MULTI-PHASE	
Recon All Images	On	Seperate Series	0
CONTRAST		Trigger Delay without AV	0
Contrast Yes/No	Yes	Mask Phase	0
		Mask Pause	0
		TRICKS	
		Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	2

+ C 3D Ax T1 FSPGR

+ C 3D Ax T1 FSPGR

Protocol: adult_head_Seizure Brain w/w/o

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	48Ch Head Coil	Imaging Options	EDR, TRF, Fast
Plane	OBLIQUE	SCAN RANGE	
Series Description	+C Ax T1 FSE	FOV	24.0
SCAN TIMING		Slice Thickness	4.0
Flip Angle	111	Slice Spacing	1.5
TE	11.0	ACQ TIMING	
Number of Echoes	1	Freq	256
TR	750.0	Phase	192
Echo Train Length	3	Freq DIR	A/P
Receiver Bandwidth	41.67	Fat Shift DIR	Normal (A)
IMAGE ENHANCE		NEX	1.00
Filter Choice	G	# of Acq. Before Pause	0
GATING/TRIGGER		Phase FOV	1.00
Auto Trigger Type	Off	Auto Shim	Auto
FMRI		Phase Correction	Yes
PSD Trigger	Internal	USER CVS	
View Order	Bottom/Up	User CV16	1.00
# of Repetitions REST	0	User CV22	1.00
# of Repetitions ACTIVE	0	TR Min	750.0
SAT		TR Max	950.0
SAT Location	I	MULTI-PHASE	
Tag Type	None	Seperate Series	0
TRICKS		Mask Phase	0
Pause On/Off	On	Mask Pause	0
Auto Subtract	0	DIFFUSION	
Auto SCIC	2	Recon All Images	On
		CONTRAST	
		Contrast Yes/No	Yes

+C Ax T1 FSE

+C Ax T1 FSE