

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
Coil Configuration	<i>AA;Body</i>	Imaging Options	<i>Seq, EDR, Fast, SS</i>
Plane	<i>3-PLANE</i>	SCAN RANGE	
Series Description	<i>3 Plane Localizer</i>	FOV	<i>48.0</i>
SCAN TIMING		Slice Thickness	<i>5.0</i>
TE	<i>80.0</i>	Slice Spacing	<i>5.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TR	<i>1204.2</i>	Freq	<i>256</i>
Receiver Bandwidth	<i>62.50</i>	Phase	<i>128</i>
IMAGE ENHANCE		Freq DIR	<i>Unswap</i>
Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>
GATING/TRIGGER		Phase FOV	<i>1.00</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
FMRI		Phase Correction	<i>No</i>
PSD Trigger	<i>Internal</i>	USER CVS	
View Order	<i>Bottom/Up</i>	User CV2	<i>240.00</i>
# of Repetitions REST	<i>0</i>	User CV13	<i>1.00</i>
# of Repetitions ACTIVE	<i>0</i>	MULTI-PHASE	
SAT		Seperate Series	<i>0</i>
Tag Type	<i>None</i>	Mask Phase	<i>0</i>
TRICKS		Mask Pause	<i>0</i>
Pause On/Off	<i>On</i>	DIFFUSION	
Auto Subtract	<i>0</i>	Recon All Images	<i>On</i>
Auto SCIC	<i>Off</i>	CONTRAST	
		Contrast Yes/No	<i>No</i>

3 Plane Localizer

3 Plane Localizer

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

Cor FIESTA ASPIR	PATIENT POSITION		IMAGING PARAMETERS		Cor FIESTA ASPIR
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>Seq, EDR, Fast, ARC</i>	
	Plane	<i>CORONAL</i>	SCAN RANGE		
	Series Description	<i>Cor FIESTA ASPIR</i>	FOV	<i>44.0</i>	
	SCAN TIMING		Slice Thickness	<i>7.0</i>	
	Flip Angle	<i>50</i>	Slice Spacing	<i>1.0</i>	
	TE	<i>Minimum</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>192</i>	
	TI	<i>210</i>	Phase	<i>300</i>	
	Receiver Bandwidth	<i>125.00</i>	Freq DIR	<i>S/I</i>	
	IMAGE ENHANCE		NEX	<i>1.00</i>	
	Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>	
	GATING/TRIGGER		Phase FOV	<i>1.00</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>	
	MULTI-PHASE		Phase Correction	<i>No</i>	
	Seperate Series	<i>0</i>	FMRI		
	Mask Phase	<i>0</i>	PSD Trigger	<i>Internal</i>	
	Mask Pause	<i>0</i>	View Order	<i>Bottom/Up</i>	
	DIFFUSION		# of Repetitions REST	<i>0</i>	
	Recon All Images	<i>On</i>	# of Repetitions ACTIVE	<i>0</i>	
	CONTRAST		SAT		
	Contrast Yes/No	<i>No</i>	Tag Type	<i>None</i>	
			Fat/Water Saturation	<i>Fat Special</i>	
		TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>Off</i>		
		OTHERS			
		Protocol Notes	<i>Kidney thru Pelvis</i>		

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, SS, Asset</i>
Plane	<i>AXIAL</i>	SCAN RANGE	
Series Description	<i>Ax SSFSE ASSET</i>	FOV	<i>42.0</i>
SCAN TIMING		Slice Thickness	<i>7.0</i>
TE	<i>90.0</i>	Slice Spacing	<i>1.0</i>
Number of Echoes	<i>1</i>	ACQ TIMING	
TR	<i>2000.0</i>	Freq	<i>320</i>
Receiver Bandwidth	<i>83.33</i>	Phase	<i>224</i>
IMAGE ENHANCE		Freq DIR	<i>R/L</i>
Filter Choice	<i>A</i>	NEX	<i>1.00</i>
GATING/TRIGGER		# of Acq. Before Pause	<i>0</i>
Auto Trigger Type	<i>Off</i>	Phase FOV	<i>1.00</i>
FMRI		Auto Shim	<i>On</i>
PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>
View Order	<i>Bottom/Up</i>	USER CVS	
# of Repetitions REST	<i>0</i>	User CV1	<i>1.00</i>
# of Repetitions ACTIVE	<i>0</i>	User CV2	<i>184.00</i>
SAT		User CV13	<i>1.00</i>
Tag Type	<i>None</i>	User CV15	<i>1.00</i>
TRICKS		MULTI-PHASE	
Pause On/Off	<i>On</i>	Seperate Series	<i>0</i>
Auto Subtract	<i>0</i>	Mask Phase	<i>0</i>
Auto SCIC	<i>2</i>	Mask Pause	<i>0</i>
OTHERS		DIFFUSION	
Protocol Notes	<i>Kidney thru Iliac crest</i>	Recon All Images	<i>On</i>
		CONTRAST	
		Contrast Yes/No	<i>No</i>

AX SSFSE ASSET

AX SSFSE ASSET

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

Ax T2 frFSE SFOV	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>FC, NPW, EDR, TRF, Fast, FR, ARC</i>
	Plane	<i>OBLIQUE</i>	SCAN RANGE	
	Series Description	<i>Ax T2 frFSE SFOV</i>	FOV	<i>22.0</i>
	SCAN TIMING		Slice Thickness	<i>3.0</i>
	Flip Angle	<i>160</i>	Slice Spacing	<i>1.0</i>
	TE	<i>102.0</i>	ACQ TIMING	
	Number of Echoes	<i>1</i>	Freq	<i>384</i>
	TR	<i>8645.0</i>	Phase	<i>256</i>
	Echo Train Length	<i>22</i>	Freq DIR	<i>A/P</i>
	Receiver Bandwidth	<i>41.67</i>	Fat Shift DIR	<i>Normal (A)</i>
	IMAGE ENHANCE		NEX	<i>1.00</i>
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
	GATING/TRIGGER		Phase Correction	<i>No</i>
	Auto Trigger Type	<i>Off</i>	Flow Direction	<i>Slice</i>
	FMRI		Compensation	
	PSD Trigger	<i>Internal</i>	USER CVS	
	View Order	<i>Bottom/Up</i>	User CV17	<i>1.00</i>
	# of Repetitions REST	<i>0</i>	User CV18	<i>1.00</i>
	# of Repetitions ACTIVE	<i>0</i>	User CV22	<i>1.00</i>
	SAT		TR Min	<i>2500.0</i>
	SAT Location	<i>R</i>	TR Max	<i>11000.0</i>
	SAT Location	<i>L</i>	MULTI-PHASE	
Tag Type	<i>None</i>	Seperate Series	<i>0</i>	
TRICKS		Mask Phase	<i>0</i>	
Pause On/Off	<i>On</i>	Mask Pause	<i>0</i>	
Auto Subtract	<i>0</i>	DIFFUSION		
Auto SCIC	<i>2</i>	Recon All Images	<i>On</i>	
CONTRAST		TRACKER		
Contrast Yes/No	<i>No</i>	Tracker Length	<i>200.0</i>	
		Tracker Thickness	<i>20.0</i>	

Ax T2 frFSE SFOV

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>FC, NPW, EDR, TRF, Fast, FR, ARC</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>Sag T2 frFSE SFOV</i>	FOV	<i>24.0</i>
SCAN TIMING		Slice Thickness	<i>3.5</i>
Flip Angle	<i>130</i>	Slice Spacing	<i>0.5</i>
TE	<i>102.0</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>384</i>
TR	<i>6452.0</i>	Phase	<i>256</i>
Echo Train Length	<i>20</i>	Freq DIR	<i>Unswap</i>
Receiver Bandwidth	<i>35.71</i>	Fat Shift DIR	<i>Normal</i>
IMAGE ENHANCE		NEX	<i>2.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
GATING/TRIGGER		Phase Correction	<i>No</i>
Auto Trigger Type	<i>Off</i>	USER CVS	
FMRI		User CV22	<i>1.00</i>
PSD Trigger	<i>Internal</i>	TR Min	<i>2200.0</i>
View Order	<i>Bottom/Up</i>	TR Max	<i>17000.0</i>
# of Repetitions REST	<i>0</i>	MULTI-PHASE	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
SAT		Mask Phase	<i>0</i>
Tag Type	<i>None</i>	Mask Pause	<i>0</i>
TRICKS		DIFFUSION	
Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>
Auto Subtract	<i>0</i>	CONTRAST	
Auto SCIC	<i>2</i>	Contrast Yes/No	<i>No</i>

Sag T2 frFSE SFOV

Sag T2 frFSE SFOV

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>FC, NPW, EDR, TRF, Fast, FR, ARC</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>Cor T2 frFSE SFOV</i>	FOV	<i>24.0</i>
SCAN TIMING		Slice Thickness	<i>3.0</i>
Flip Angle	<i>130</i>	Slice Spacing	<i>1.0</i>
TE	<i>102.0</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>384</i>
TR	<i>6410.0</i>	Phase	<i>256</i>
Echo Train Length	<i>20</i>	Freq DIR	<i>Unswap</i>
Receiver Bandwidth	<i>35.71</i>	Fat Shift DIR	<i>Normal</i>
IMAGE ENHANCE		NEX	<i>2.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
GATING/TRIGGER		Phase Correction	<i>No</i>
Auto Trigger Type	<i>Off</i>	USER CVS	
FMRI		User CV22	<i>1.00</i>
PSD Trigger	<i>Internal</i>	TR Min	<i>2200.0</i>
View Order	<i>Bottom/Up</i>	TR Max	<i>17000.0</i>
# of Repetitions REST	<i>0</i>	MULTI-PHASE	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
SAT		Mask Phase	<i>0</i>
Tag Type	<i>None</i>	Mask Pause	<i>0</i>
TRICKS		DIFFUSION	
Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>
Auto Subtract	<i>0</i>	CONTRAST	
Auto SCIC	<i>2</i>	Contrast Yes/No	<i>No</i>

Cor T2 frFSE SFOV

Cor T2 frFSE SFOV

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

Ax DWI All B-50/500/800	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Spin Echo</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EPI, DIFF, Asset</i>
	Plane	<i>AXIAL</i>	SCAN RANGE	
	Series Description	<i>Ax DWI All B-50/500/800</i>	FOV	<i>26.0</i>
	SCAN TIMING		Slice Thickness	<i>4.0</i>
	TE	<i>Minimum</i>	Slice Spacing	<i>0.0</i>
	Number of Echoes	<i>1</i>	ACQ TIMING	
	TR	<i>6588.0</i>	Freq	<i>96</i>
	Number of Shots	<i>1</i>	Phase	<i>96</i>
	IMAGE ENHANCE		Freq DIR	<i>R/L</i>
	Filter Choice	<i>None</i>	Phase FOV	<i>1.00</i>
	GATING/TRIGGER		Auto Shim	<i>Auto</i>
	Auto Trigger Type	<i>Off</i>	Phase Correction	<i>Yes</i>
	FMRI		USER CVS	
	PSD Trigger	<i>Internal</i>	User CV0	<i>1.00</i>
	View Order	<i>Bottom/Up</i>	User CV5	<i>1.00</i>
	# of Repetitions REST	<i>0</i>	User CV17	<i>1.00</i>
	# of Repetitions ACTIVE	<i>0</i>	User CV18	<i>1.00</i>
	SAT		TR Min	<i>3000.0</i>
	SAT Location	<i>A</i>	TR Max	<i>17000.0</i>
	Tag Type	<i>None</i>	MULTI-PHASE	
	Fat/Water Saturation	<i>Fat</i>	Seperate Series	<i>0</i>
	TRICKS		Mask Phase	<i>0</i>
Pause On/Off	<i>On</i>	Mask Pause	<i>0</i>	
Auto Subtract	<i>0</i>	DIFFUSION		
Auto SCIC	<i>2</i>	Optimized TE	<i>Yes</i>	
		Diffusion Directions	<i>All</i>	
		Number of Diffusion Directions	<i>3</i>	
		Number of T2 Images	<i>0</i>	
		Dual Spin Echo	<i>Off</i>	
		Recon All Images	<i>On</i>	
		CONTRAST		
		Contrast Yes/No	<i>No</i>	

Ax DWI All B-50/500/800

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>NPW, EDR, Fast, ARC</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>Ax T1 FSPGR</i>	FOV	<i>35.0</i>
SCAN TIMING		Slice Thickness	<i>4.0</i>
Flip Angle	<i>60</i>	Slice Spacing	<i>0.5</i>
TE	<i>In Phase</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>384</i>
TR	<i>450.0</i>	Phase	<i>224</i>
Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>R/L</i>
IMAGE ENHANCE		NEX	<i>2.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
GATING/TRIGGER		Phase Correction	<i>No</i>
Auto Trigger Type	<i>Off</i>	FMRI	
MULTI-PHASE		PSD Trigger	<i>Internal</i>
Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
DIFFUSION		SAT	
Recon All Images	<i>On</i>	SAT Location	<i>S</i>
CONTRAST		SAT Location	<i>I</i>
Contrast Yes/No	<i>No</i>	Tag Type	<i>None</i>
		TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>2</i>

Ax T1 FSPGR

Ax T1 FSPGR

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>NPW, EDR, Fast, ARC</i>
Plane	<i>OBLIQUE</i>	SCAN RANGE	
Series Description	<i>Ax T1 FSPGR F/S</i>	FOV	<i>35.0</i>
SCAN TIMING		Slice Thickness	<i>4.0</i>
Flip Angle	<i>60</i>	Slice Spacing	<i>0.5</i>
TE	<i>In Phase</i>	ACQ TIMING	
Number of Echoes	<i>1</i>	Freq	<i>384</i>
TR	<i>500.0</i>	Phase	<i>224</i>
Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>A/P</i>
IMAGE ENHANCE		NEX	<i>2.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
GATING/TRIGGER		Phase Correction	<i>No</i>
Auto Trigger Type	<i>Off</i>	FMRI	
MULTI-PHASE		PSD Trigger	<i>Internal</i>
Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>
Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>
Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
DIFFUSION		SAT	
Recon All Images	<i>On</i>	SAT Location	<i>S</i>
CONTRAST		SAT Location	<i>I</i>
Contrast Yes/No	<i>No</i>	Tag Type	<i>None</i>
		Fat/Water Saturation	<i>Fat</i>
		TRICKS	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>2</i>

Ax T1 FSPGR F/S

Ax T1 FSPGR F/S

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	Cube
Coil Configuration	Body	Imaging Options	EDR, Fast, ZIP512, ZIP2, FR, ARC
Plane	AXIAL		
Series Description	3D Ax T2 HyperCube		
SCAN TIMING		SCAN RANGE	
TE	99.0	FOV	30.0
Number of Echoes	1	Slice Thickness	2.4
TR	2000.0	Location per Slab	76
Echo Train Length	70	Overlap Locations	0
Receiver Bandwidth	41.67		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	256
		Phase	256
		Freq DIR	R/L
		NEX	1.00
		Phase FOV	0.80
		Auto Shim	Auto
		Phase Correction	No
GATING/TRIGGER		FMRI	
Auto Trigger Type	Off	PSD Trigger	Internal
		View Order	Bottom/Up
		# of Repetitions REST	0
		# of Repetitions ACTIVE	0
MULTI-PHASE		SAT	
Seperate Series	0	Tag Type	None
Mask Phase	0		
Mask Pause	0		
DIFFUSION		TRICKS	
Recon All Images	On	Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	2
CONTRAST			
Contrast Yes/No	No		

3D Ax T2 HyperCube

3D Ax T2 HyperCube

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, Asset</i>
Plane	<i>AXIAL</i>	SCAN RANGE	
Series Description	<i>3D Ax LAVA Pre</i>	FOV	<i>32.0</i>
SCAN TIMING		Slice Thickness	<i>4.5</i>
Flip Angle	<i>10</i>	Location per Slab	<i>32</i>
Number of Echoes	<i>1</i>	Overlap Locations	<i>0</i>
TI	<i>5</i>	ACQ TIMING	
Receiver Bandwidth	<i>62.50</i>	Freq	<i>192</i>
IMAGE ENHANCE		Phase	<i>160</i>
Filter Choice	<i>None</i>	Freq DIR	<i>A/P</i>
GATING/TRIGGER		Phase FOV	<i>0.90</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
FMRI		Phase Correction	<i>No</i>
PSD Trigger	<i>Internal</i>	USER CVS	
View Order	<i>Bottom/Up</i>	User CV6	<i>1.00</i>
# of Repetitions REST	<i>0</i>	MULTI-PHASE	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
SAT		Delay after Acquisition without AV	<i>1</i>
Tag Type	<i>None</i>	Trigger Delay without AV	<i>0</i>
Fat/Water Saturation	<i>Fat Special</i>	Mask Phase	<i>0</i>
TRICKS		Mask Pause	<i>0</i>
Pause On/Off	<i>On</i>	DIFFUSION	
Auto Subtract	<i>0</i>	Recon All Images	<i>On</i>
Auto SCIC	<i>2</i>	CONTRAST	
		Contrast Yes/No	<i>No</i>

3D Ax LAVA Pre

3D Ax LAVA Pre

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA</i>
Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, Asset</i>
Plane	<i>AXIAL</i>	SCAN RANGE	
Series Description	<i>+c 3D Ax LAVA</i>	FOV	<i>32.0</i>
SCAN TIMING		Slice Thickness	<i>4.5</i>
Flip Angle	<i>10</i>	Location per Slab	<i>32</i>
Number of Echoes	<i>1</i>	Overlap Locations	<i>0</i>
TI	<i>5</i>	ACQ TIMING	
Receiver Bandwidth	<i>62.50</i>	Freq	<i>192</i>
IMAGE ENHANCE		Phase	<i>160</i>
Filter Choice	<i>None</i>	Freq DIR	<i>A/P</i>
GATING/TRIGGER		Phase FOV	<i>0.90</i>
Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>
FMRI		Phase Correction	<i>No</i>
PSD Trigger	<i>Internal</i>	USER CVS	
View Order	<i>Bottom/Up</i>	User CV6	<i>1.00</i>
# of Repetitions REST	<i>0</i>	MULTI-PHASE	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
SAT		Delay after Acquisition without AV	<i>1</i>
Tag Type	<i>None</i>	Trigger Delay without AV	<i>0</i>
Fat/Water Saturation	<i>Fat Special</i>	Mask Phase	<i>0</i>
TRICKS		Mask Pause	<i>0</i>
Pause On/Off	<i>On</i>	DIFFUSION	
Auto Subtract	<i>0</i>	Recon All Images	<i>On</i>
Auto SCIC	<i>2</i>	CONTRAST	
		Contrast Yes/No	<i>Yes</i>

+c 3D Ax LAVA

+c 3D Ax LAVA

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

+c Ax T1 FSPGR F/S SFOV	PATIENT POSITION		IMAGING PARAMETERS		+c Ax T1 FSPGR F/S SFOV
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>NPW, EDR, Fast, ARC</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>+c Ax T1 FSPGR F/S SFOV</i>	FOV	<i>22.0</i>	
	SCAN TIMING		Slice Thickness	<i>3.0</i>	
	Flip Angle	<i>60</i>	Slice Spacing	<i>1.0</i>	
	TE	<i>In Phase</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>384</i>	
	TR	<i>500.0</i>	Phase	<i>224</i>	
	Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>Unswap</i>	
	IMAGE ENHANCE		NEX	<i>2.00</i>	
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>	
	GATING/TRIGGER		Phase Correction	<i>No</i>	
	Auto Trigger Type	<i>Off</i>	FMRI		
	MULTI-PHASE		PSD Trigger	<i>Internal</i>	
	Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>	
	Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>	
	Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>	
	DIFFUSION		SAT		
	Recon All Images	<i>On</i>	SAT Location	<i>S</i>	
	CONTRAST		SAT Location	<i>I</i>	
	Contrast Yes/No	<i>Yes</i>	Tag Type	<i>None</i>	
			Fat/Water Saturation	<i>Fat</i>	
		TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>2</i>		

McLeod Pavilion

Protocol: adult_pelvis_Conwell Female Pelvis

+c Sag T1 FSPGR F/S SFOV	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
	Coil Configuration	<i>Body</i>	Imaging Options	<i>NPW, EDR, Fast, ARC</i>
	Plane	<i>OBLIQUE</i>	SCAN RANGE	
	Series Description	<i>+c Sag T1 FSPGR F/S SFOV</i>	FOV	<i>24.0</i>
	SCAN TIMING		Slice Thickness	<i>3.5</i>
	Flip Angle	<i>60</i>	Slice Spacing	<i>0.5</i>
	TE	<i>In Phase</i>	ACQ TIMING	
	Number of Echoes	<i>1</i>	Freq	<i>320</i>
	TR	<i>375.0</i>	Phase	<i>224</i>
	Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>Unswap</i>
	IMAGE ENHANCE		NEX	<i>2.00</i>
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
	GATING/TRIGGER		Phase Correction	<i>No</i>
	Auto Trigger Type	<i>Off</i>	FMRI	
	MULTI-PHASE		PSD Trigger	<i>Internal</i>
	Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>
	Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>
	Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
	DIFFUSION		SAT	
	Recon All Images	<i>On</i>	Tag Type	<i>None</i>
	CONTRAST		Fat/Water Saturation	<i>Fat</i>
	Contrast Yes/No	<i>Yes</i>	TRICKS	
		Pause On/Off	<i>On</i>	
		Auto Subtract	<i>0</i>	
		Auto SCIC	<i>2</i>	

+c Sag T1 FSPGR F/S SFOV

Protocol: adult_pelvis_Conwell Female Pelvis

+c Cor T1 FSPGR F/S SFOV	PATIENT POSITION		IMAGING PARAMETERS		+c Cor T1 FSPGR F/S SFOV
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>NPW, EDR, Fast, ARC</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>+c Cor T1 FSPGR F/S SFOV</i>	FOV	<i>24.0</i>	
	SCAN TIMING		Slice Thickness	<i>3.0</i>	
	Flip Angle	<i>60</i>	Slice Spacing	<i>1.0</i>	
	TE	<i>In Phase</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>320</i>	
	TR	<i>450.0</i>	Phase	<i>224</i>	
	Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>S/I</i>	
	IMAGE ENHANCE		NEX	<i>2.00</i>	
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>	
	GATING/TRIGGER		Phase Correction	<i>No</i>	
	Auto Trigger Type	<i>Off</i>	FMRI		
	MULTI-PHASE		PSD Trigger	<i>Internal</i>	
	Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>	
	Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>	
	Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>	
DIFFUSION		SAT			
Recon All Images	<i>On</i>	Tag Type	<i>None</i>		
CONTRAST		Fat/Water Saturation	<i>Fat</i>		
Contrast Yes/No	<i>Yes</i>	TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>2</i>		

Protocol: adult_pelvis_Conwell Female Pelvis

+c Ax T1 FSPGR F/S LFOV	PATIENT POSITION		IMAGING PARAMETERS		+c Ax T1 FSPGR F/S LFOV
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>	
	Coil Configuration	<i>Body</i>	Imaging Options	<i>EDR, Fast, ARC</i>	
	Plane	<i>OBLIQUE</i>	SCAN RANGE		
	Series Description	<i>+c Ax T1 FSPGR F/S LFOV</i>	FOV	<i>40.0</i>	
	SCAN TIMING		Slice Thickness	<i>5.0</i>	
	Flip Angle	<i>60</i>	Slice Spacing	<i>1.0</i>	
	TE	<i>In Phase</i>	ACQ TIMING		
	Number of Echoes	<i>1</i>	Freq	<i>384</i>	
	TR	<i>400.0</i>	Phase	<i>192</i>	
	Receiver Bandwidth	<i>31.25</i>	Freq DIR	<i>R/L</i>	
	IMAGE ENHANCE		NEX	<i>2.00</i>	
	Filter Choice	<i>None</i>	Phase FOV	<i>1.00</i>	
	GATING/TRIGGER		Auto Shim	<i>Auto</i>	
	Auto Trigger Type	<i>Off</i>	Phase Correction	<i>No</i>	
	MULTI-PHASE		FMRI		
	Seperate Series	<i>0</i>	PSD Trigger	<i>Internal</i>	
	Mask Phase	<i>0</i>	View Order	<i>Bottom/Up</i>	
	Mask Pause	<i>0</i>	# of Repetitions REST	<i>0</i>	
	DIFFUSION		# of Repetitions ACTIVE	<i>0</i>	
	Recon All Images	<i>On</i>	SAT		
	CONTRAST		SAT Location	<i>S</i>	
	Contrast Yes/No	<i>No</i>	SAT Location	<i>I</i>	
			Tag Type	<i>None</i>	
		Fat/Water Saturation	<i>Fat</i>		
		TRICKS			
		Pause On/Off	<i>On</i>		
		Auto Subtract	<i>0</i>		
		Auto SCIC	<i>2</i>		